

**Lab [25 pts]**  
**Molarity and Absorption**

Name \_\_\_\_\_  
Lab Partner \_\_\_\_\_  
Period \_\_\_\_\_ Date \_\_\_\_\_

**Purpose:** To become familiar with making solutions by dissolving solids into water and by diluting solutions. To experimentally verify that the concentration of a solution is directly proportional to the amount of light absorbed (absorbance) in a colorimeter (Beer's Law).

**Background:** Chemists often need to determine the concentration of a solute in a solution. They can do this by using a spectrophotometer or colorimeter to measure the amount of light of a given wavelength absorbed by the solution. By choosing a wavelength of light of the complimentary color to the solution, the absorbance of the solution is directly proportional to the concentration, a relationship called **Beer's Law**:  $A = a \times b \times c$  where A is the absorbance (what you measure), a is the absorptivity of the solute, b is the path length of the cuvette (see below), and c is the concentration. For us, a and b will be constant, so  $A = mc$ , where m is the slope of the line formed by plotting A for set of solutions of known concentration, called a calibration "curve" (even though it should be a straight line). From this, you can determine an unknown concentration by measuring A.

**Prelab Calculations/Procedure:** Making solutions of  $\text{CuSO}_4$  and testing their absorbance. *Perform all indicated calculations before coming into the lab!*

- 1) **0.50 M Solution: Make 50 mL of 0.50 M  $\text{CuSO}_4$  solution from solid  $\text{CuSO}_4$  crystals**
  - a) [2 pt Prelab] Determine the grams of  $\text{CuSO}_4$  needed to make 50 mL of 0.50 M  $\text{CuSO}_4$  (aq).  
\*\* **Needed info:** *The molar mass is 249.6 g/mole since  $\text{CuSO}_4(s)$  is a hydrate, which has 5 water molecules contained within the crystal lattice of the solid. Its formula is  $\text{CuSO}_4 \cdot 5 \text{H}_2\text{O}$ . Show all work here!*
  - b) Place the needed mass of  $\text{CuSO}_4$  (from part a) into a **50 mL volumetric flask** (Careful, don't spill). Fill about 3/4 with distilled  $\text{H}_2\text{O}$ . Shake flask until all of the solid dissolves. Be patient!!! Add distilled water up to the 50 mL mark, cap flask and shake until well mixed. Pour solution into beaker.
  - c) Fill a plastic cuvette with the 0.50 M solution you just made. Measure its absorbance and record in data table on page 2, below its concentration. *See directions on next page to learn how to measure absorbance. Record in Data Table on next page.*
- 2) **0.35 M Solution: Make 25 mL of a 0.35 M  $\text{CuSO}_4$  solution by diluting your 0.5 M solution from #1.**
  - a) [2 pt Prelab] Use the dilution equation to determine the volume of the 0.5 M  $\text{CuSO}_4$  solution that must be diluted with water to make 25 mL of 0.35 M  $\text{CuSO}_4$ . Show all work here!
  - b) Use a **10-mL graduated pipette** to measure out the needed volume (from part a) of your 0.50 M  $\text{CuSO}_4$ . Pour solution into a **25 mL volumetric flask**. Add distilled water up to the 25 mL mark. Cap the flask and shake until well mixed.
  - c) Fill a plastic cuvette with the 0.35 M solution you just made. Measure its absorbance and record in data table on page 2, below its concentration.
- 3) **Making four solutions by dilution of your 0.50 M  $\text{CuSO}_4$  solution:** a simpler method to make dilute solutions.
  - a) [2 pt Prelab] Use the dilution equation to calculate the volume of 0.50 M solution (which is  $M_1$  for all) needed to make 10.0 mL ( $V_2$ ) of the 4 dilutions as shown in the table below. To find the volume of water needed, subtract  $V_1$  from 10.0 mL. **Show sample calculation for test tube #1 below table.**

Test tube	#1	#2	#3	#4
Molarity of solution, M ( $M_2$ )	0.10	0.20	0.30	0.40
Vol of 0.5 M $\text{CuSO}_4$ , mL ( $V_1$ )				
Vol of water needed, mL				

Sample calculations for Test Tube #1:

- Label four large test tubes #1, #2, #3, and #4. Place them in your test tube rack. Use two **10-mL graduated pipettes** to make the 4 different dilutions by filling each test tube with the amounts of each solution shown above. Put the  $\text{CuSO}_4$  solution in each test tube first, then add the water.
- Measure the absorbance of each of the solutions just made. *Record in Data Table, below proper concentration.*

\*\*\* **To measure Absorbance:** Use the colorimeter/computer setup:

- Carefully fill a colorimeter **cuvette** (plastic container) with one of the solutions. **Cap and wipe dry!!!**
- Take the filled cuvette to the colorimeter/computer. Put the sample into the sample chamber. (*Make sure that the clear sides face the white line on the colorimeter.*) Close the lid of the colorimeter. Read the **absorbance** off the computer screen. Record value in data table below concentration of solution.
- Remove cuvette. Go back to lab bench. Pour solution down drain, rinse with a new solution, refill with the same new solution, and repeat absorbance measurement.
- When you are done, empty 0.50 M solution into beaker at front, all other solutions can be rinsed down the sink, and all of your equipment *must* be rinsed with deionized  $\text{H}_2\text{O}$  until no blue remains.

**Data Table** [2 pts]

Molarity (M)	0.10 M	0.20 M	0.30 M	0.35 M	0.40 M	0.50 M
Absorbance						

**Graph:** [6 pts] **Beer's Law Plot of Absorbance vs. Molarity**

- Make a graph of the absorbance vs. molarity of your six solutions. Remember-- labels, units, and title. *Absorbance goes on the y-axis and Molarity goes on the x-axis.*
- You should have six points on your graph and they should lie in a straight line. **Draw the best fit straight line** (use a ruler!) through the six points, making sure that it goes through the origin (0,0) since a solution of 0 M should have  $A = 0$ !
- Make sure you show all labels, units and have a descriptive title (Not just absorbance vs. Molarity)

**Post-Lab Questions:** [11 pts] *Answer on separate sheet.*

- [2 pts] Is absorbance directly or inversely related to molarity? How do you know?
- [1 pt] Suppose you have a solution of  $\text{CuSO}_4$  of unknown molarity. You measure its absorbance and it turns out to be 0.55. What is the molarity of the solution based on your absorbance data obtained in lab?  
\*\* *Look at your best-fit line and determine the molarity on the line that corresponds to a 0.55 absorbance.*
  - If the absorbance is 0.55, the molarity = \_\_\_\_\_
  - On your graph, show where the point is by drawing an asterisk (\*) at the correct spot on the line.
- [2 pts] By looking at your data, you know that absorption of light increases with increasing concentration. How would you have known this just by looking at the actual solutions that you made? [How would you describe the color of the 0.5 M solution compared to the 0.1 M solution?]
- [3 pts] Suppose you wanted to make 100.0 mL of a 0.500 M  $\text{CaCl}_2$  solution from solid  $\text{CaCl}_2$ .
  - Calculate how many grams of  $\text{CaCl}_2$  will be needed. Show **all steps and units** and **watch sig figs!!**
  - Briefly describe how you would make this solution in the lab. (*Be specific about amounts of substances and the size and type of glassware!*)
- [3 pts] Suppose you wish to make a 225 mL of 2.00 M  $\text{HNO}_3$  (aq). You will accomplish your goal by diluting 15.0 M  $\text{HNO}_3$  (aq).
  - What volume of the 15.0 M  $\text{HNO}_3$  (aq) solution is needed? *Watch units and sig figs!!*
  - Briefly describe how you would make this solution in the lab. (*Be specific about amounts of substances and the size and type of glassware!*)